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Principle : Methylene blue is a simple stain that is particularly useful in the identification of *Corynebacterium* species.

The metachromatic granules of *C. diphtheriae* readily take up Methylene blue dye and appear deep blue. Although some authors have stated that the cytoplasmic granule formation characteristic of *C. diphtheriae* is rarely seen with saprophytic species of *Corynebacterium*, this criterion is unreliable and cannot be used for definitive identification of *C. diphtheriae* without further studies.

Procedure

1. Heat-fix the smear. Flood the surface of the smear with the Methylene blue-staining solution for 1 minute. Wash the slide with water and blot dry.
2. In the past it was necessary to add alkali to the above solution before use. However, Methylene blue dyes prepared in recent years do not require this additional step because acid impurities found in older stains have been removed.

Results

Interpretation

The corynebacteria are pleomorphic bacilli that range in size from 0.5 to 1.0 μm in width and from 2 to 6 μm in length and appear as straight, curved, or club-shaped rods. Characteristic for the micro-organisms are metachromatic granules that take up Methylene blue stain and appear dark blue. Although this finding is characteristic of the corynebacteria, species of the

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TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI\RESP\11\v05	Page 1 of 5
Section: Respiratory Tract Culture Manual	Subject Title: SPUTUM (Including Endotracheal Tube and Tracheostomy Specimens)	
Issued by: LABORATORY MANAGER	Original Date: September 25, 2000	
Approved by: Laboratory Director	Revision Date: September 14, 2006	
	Annual Review Date: August 13, 2007	

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TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI\RESP\11\v05	Page 2 of 5
Respiratory Tract Culture Manual		

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CULTURES

- URINES** — 1/2 BA & CLED or MAC + 1/2 BA
Clinical data Typhoid, Pso & gastro-enteritis — Centrifuge & put day into Selenite F broth overnight then culture.
- STOOL** — FDCA or XLD, Selenite F broth — incubate overnight.
Under 2 yrs — EMB, BA or MAC
ALL rice water stools — TCBS & alkaline peptone water for vibrios
- BLOOD** — incubate at 37°C overnight.
Subculture 10 days — 1st, 5th, 10th day. BA (anO₂), MAC & CHOC
- VAGINAL SWAB (cervical)** — Wet prep: epithelial cells, wbc's, T. va. clue cells & spermatozoa.
Gram: clue cells and organisms
BA, CHOC & MAC.
(TM or NYC) for gonococcal isolation
Sweep choc.
as vaginal swab.
- URETHRAL/PENILE DISCHARGE** — as vaginal swab.
- THROAT SWAB** — Gram stain - BA, CHOC & MAC
- EAR SWAB** — Wet Prep & Gram Stain. BA, CHOC, MAC
- EYE SWAB** — Same as ear.
- WOUNDS & Fluids** — All surgical specimens eg Catheter tips & in Robertson's medium, incubate at 37°C overnight
BA, CHOC, MAC
- Burns**
- Fluids** — wet prep. 1/2 BA, CHOC, MAC & AZIDS
Deposit — Gram. Leish, ZN & in Robertson
1/2 BA, CHOC, MAC. All joint fluids for culture.

AFB SMEAR STAINING

- Always use new, grease free, and clean slides. Correctly label slides with stylus or lead pencil.
- Fish out yellowish portion from sputum container and place on slide with the rough end of the stick.
- Spread material evenly in an approximate area of 2cm X 1cm so that new point is visible on drying.
- Air dry smear completely and then heat fix smear in a flame.
- Place slides on the staining rack without touching each other. Always add Positive and Negative control slides.
- Cover slides with freshly filtered carbol fuchsin.
- Heat gently with a torch until steam rises from the slides. Stain for five minutes.
- Wash gently with water.
- Drain the water.
- Cover slides with decolorizing solution for three minutes.
- Wash thoroughly with water. If slide is not decolorized properly repeat step 10 for additional 1-3 minutes. Rinse thoroughly with water.
- Drain the water.
- Cover with counter stain. Methylene blue for one minute.
- Drain the counter stain.
- Wash with water. Wipe the back side of slides with tissue paper.
- Air dry the slides in a rack.
- View the smear under oil immersion. AFB: Fine, red rods against blue background.

AFB Counts	Recording/Reporting
0 to 9 AFB in 100 fields	—
1 to 9 AFB in 100 fields	Actual AFB count
10 to 99 AFB in 100 fields	—
1 to 10 AFB per field	—
20 to 100 AFB per field	—
> 10 AFB per field	—
at least 20 fields	—

Report the findings as per WHO and IGHUO recommendations.

A joint effort of:

The illustrations 1 and 2-17 are used with the permission of RIT/IGUO from "TB Bacteriology Examination to Step 18 by Akiko Fujita"

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